

BJ TRADE+LINK



PRODUCT CATALOGUE 2023

TABLE OF CONTENTS

Splints 2



SPLINTS



SPLINT STRAIGHT (WOODEN)

Straight wooden leg splint.


Wooden Splints are used for immobilising fractured joints and bones. Used together with a plastic junction.

Order 2 wooden splints and 1 plastic junction to make up a full leg splint.

[Read More](#)

SKU: SPL014

Price:



SPLINT PLASTIC JUNCTION

Plastic Junction for wooden splint.

Plastic junctions are used to connect 2 wooden splints. Wooden Splints are used for immobilising fractured joints and bones.

[Read More](#)

SKU: SPL011

Price:



SPLINT KRAMER (1M)

Kramer Wire Splint

The Kramer Wire Splint is designed to immobilise the injured body part and protects it from further injury. It aids in the reduction of pain, discomfort and swelling by limiting the movement of the injured limb. The Kramer Wire

Splint is an extremely pliable multi-purpose splint and is built from a thin core of aluminium wire mesh, sandwiched between two layers of closed...

[Read More](#)

SKU: SP006

Price:



SPLINT KRAMER (0.5M)

Kramer Wire Splint



The Kramer Wire Splint is designed to immobilise the injured body part and protects it from further injury. It aids in the reduction of pain, discomfort and swelling by limiting the movement of the injured limb. The Kramer Wire Splint is an extremely pliable multi-purpose splint and is built from a thin core of aluminium wire mesh, sandwiched between two layers of closed...

[Read More](#)

SKU: SP005

Price:



SPLINT INTERCHANGEABLE INTERLOCKING LEG/ARM YELLOW

Interchangeable Leg/Arm and Interlocking Yellow Splint.

Designed to protect broken bones, injured tendons, ligaments and dislocated joints by limiting movement on the injured area. This helps reduce pain, swelling and further damage.

Locks into a 180° and 90° configuration
1 tongue 2 grooves locking system

Multiple units can be added to get the desired length by clipping the tongue...

[Read More](#)

SKU: SPL004

Price:



SPLINT INTERCHANGEABLE INTERLOCKING LEG/ARM BLACK

Interchangeable Leg/Arm and Interlocking Black Splint.

Designed to protect broken bones, injured tendons, ligaments and dislocated joints by limiting movement on the injured area. This helps reduce pain, swelling and further damage.

Locks into a 180° and 90° configuration

1 tongue 2 grooves locking system

Multiple units can be added to get the desired length by clipping the tongue...

[Read More](#)

SKU: SPL003

Price:



SPLINT ELBOW (WOODEN)

Wooden elbow splint in L-shape.

Lightweight, durable, smooth wooden splints for immobilization of arms, ankles and feet.

Provides firm support for IV procedures. Can be used with splint padding or abdominal pads for comfort

Made of durable, lightweight hardboard.

[Read More](#)**SKU:** SPL001**Price:**

SPLINT COMBO ARM/LEG DISPOSABLE



Combo arm/leg 90-degree disposable splint.

An easily disposable and lightweight splint made of recyclable plastic. The combo splint can clip into a 180° and 90° configuration. Cost-effective solution. Can be used in first aid kits.

[Read More](#)**SKU:** SPL002**Price:**

SPLINT ALUMINIUM FLEXIBLE FOLD-UP



Folded Aluminium splint.

The aluminium splint bends easily allowing for various applications. The two layers of closed-cell foam and its aluminium interior can be cut with ordinary scissors to adapt to any size requirements. Can be easily cleaned and disinfected with bleach and water for reuse. Compatible with X-rays and CT-scans without removing the splint from the fractured or injured...

[Read More](#)**SKU:** SPL017



INDEX

S

- Splint Aluminium Flexible Fold-Up 6
- Splint Combo arm/leg disposable 6
- Splint Elbow (Wooden) 5
- Splint Interchangeable Interlocking Leg/Arm Black 5
- Splint Interchangeable Interlocking Leg/Arm Yellow 4
- Splint Kramer (0.5m) 4
- Splint Kramer (1m) 3
- Splint Plastic Junction 3
- Splint Straight (Wooden) 3