

# BJ TRADE+LINK



## PRODUCT CATALOGUE 2023

# TABLE OF CONTENTS

<b>First Aid Kit Items</b> .....	2
<b>First Aid Kits</b> .....	5



## **FIRST AID KIT ITEMS**

## RESUSAID CPR RED WITH CLEAR APRON



Resusaid Red CPR 165x165mm with clear apron.

A Resusaid or CPR mouthpiece, is a small device used to safely deliver rescue breaths during a cardiac arrest or respiratory arrest. Air is administered to the patient when the emergency responder exhales through the one-way filter valve. The CPR Mouthpiece has a built in one-way valve with a filter to protect the emergency responder from the...

[Read More](#)

**SKU:** RESC034

**Price:**



## RESUSAID CPR 165X165MM WITH MILKY APRON



Resusaid CPR 165x165mm with milky apron.

A Resusaid or CPR mouthpiece, is a small device used to safely deliver rescue breaths during a cardiac arrest or respiratory arrest. Air is administered to the patient when the emergency responder exhales through the one-way filter valve. The CPR Mouthpiece has a built in one-way valve with a filter to protect the emergency responder from the...

[Read More](#)

**SKU:** RESC006

**Price:**

## RESUSAID CPR 165X165MM WITH CLEAR APRON



Resusaid CPR 165x165mm with clear apron.

A Resusaid or CPR mouthpiece, is a small device used to safely deliver rescue breaths during a cardiac arrest or respiratory arrest. Air is administered to the patient when the emergency responder exhales through the one-way filter valve. The CPR Mouthpiece has a built in one-way valve with a filter to protect the emergency responder from the...

[Read More](#)

**SKU:** RESC023

**Price:**



## FIRST AID KITS

Qty Description 1 x Resusaid 165x165mm Milky with instructions 1 x Gloves Nitrile Powder-Free (Large) pair 1 x Antiseptic Wipes (200mm x 190mm) [Read More](#)

**SKU:** BAGFA006

**Price:**



# INDEX

## **R**

Resusaid CPR 165x165mm with clear apron 4

Resusaid CPR 165x165mm with Milky apron 3

Resusaid CPR Red with Clear apron 3