

# BJ TRADE+LINK



## PRODUCT CATALOGUE 2023

# TABLE OF CONTENTS

**Tourniquets** ..... 2



## **TOURNIQUETS**

## TOURNIQUET - ELASTIC & VELCRO MEDIUM BLACK



Tourniquet - Elastic & Velcro Medium Black

The tourniquet allows for pressure to be applied to the arm so that venous blood returning to the heart can be slowed down. As a result, the blood vessel walls become temporarily occluded and the veins distend due to the pooling of blood. This allows veins to become more visible and easier to palpate.

[Read More](#)

**SKU:** TOURN002

**Price:**



## TOURNIQUET - ELASTIC & BUCKLE LARGE BLACK



Tourniquet - Elastic & Buckle Large Black 45cm.

The tourniquet allows for pressure to be applied to the arm so that venous blood returning to the heart can be slowed down. As a result, the blood vessel walls become temporarily occluded and the veins distend due to the pooling of blood. This allows veins to become more visible and easier to palpate.

[Read More](#)

**SKU:** TOURN003

**Price:**





# INDEX

## T

Tourniquet - Elastic & Buckle LARGE Black 3

Tourniquet - Elastic & Velcro Medium Black 3