

BJ TRADE+LINK



PRODUCT CATALOGUE 2023

TABLE OF CONTENTS

Spine boards / Stretcher Accessories 2



SPINE BOARDS / STRETCHER ACCESSORIES

HEAD IMMOBILISER BOYD



Head Immobiliser Boyd

Easy to use and compatible with most spine boards on the market. The side blocks with wide Velcro closures, allow for a comfortable fit. It also comes with a head/chin straps and a base plate that fastens easily to the board. The Head Immobilizer has no metal components and can be worn during x-rays, MRI and CT scan procedures. Easy to clean and disinfect for...

Head Immobiliser Boyd. Easy to use and compatible with most spine boards on the market. The side blocks with wide Velcro closures, allow for a comfortable fit. It also comes with a head/chin straps and a base plate that fastens easily to the board. The Head Immobilizer has no metal components and can be worn during x-rays, MRI and CT scan procedures. Easy to clean and disinfect for reuse. The Head Immobilizer is made with the highest grade of foam.

[Read More](#)

SKU: RESC010

Price:

INDEX

H

Head Immobiliser Boyd 3