

BJ TRADE+LINK



PRODUCT CATALOGUE 2023

TABLE OF CONTENTS

First Aid Bags 2

First Aid Kit Items 4

Gloves 6

Tourniquets 9

Trauma Collars 11



FIRST AID BAGS

BAG ONLY - BLUE TECHNICIAN BACKPACK 45X36X15CM (ESKOM)



Bag only - Blue Eskom backpack for outdoors and technicians.

Dimensions: 45x36x15cm

[Read More](#)


SKU: BAGFA026

Price:





FIRST AID KIT ITEMS



SAFETY PINS LARGE SILVER 45MM - 12S IN BUNCH

12 Silver Large #3 (45mm) Safety pins in the bunch.

All same sizes.

[Read More](#)

SKU: SPIN006

Price:



GLOVES



GLOVES LATEX POWDERED LARGE (PAIRS)

Gloves Latex Powdered Large (Pairs)

[Read More](#)

SKU: LGNS009

Price:



GLOVES LATEX POWDERED LARGE (BOX 100)

Gloves Latex Powdered Large (Box 100)

[Read More](#)

SKU: LGNS001

Price:



GLOVES LATEX PF LARGE (BOX 100)

Gloves Latex Powder Free Large (Box 100)

[Read More](#)

SKU: LGNSPF003

Price:



GLOVES LATEX HIGH RISK SAFETY PF LARGE (50 PAIRS)

Gloves Latex High Risk Safety Powder-Free Large (50 pairs)

[Read More](#)

SKU: LGNS005

Price:





TOURNIQUETS

TOURNIQUET - ELASTIC & BUCKLE LARGE BLACK



Tourniquet - Elastic & Buckle Large Black 45cm.

The tourniquet allows for pressure to be applied to the arm so that venous blood returning to the heart can be slowed down. As a result, the blood vessel walls become temporarily occluded and the veins distend due to the pooling of blood. This allows veins to become more visible and easier to palpate.

[Read More](#)

SKU: TOURN003

Price:





TRAUMA COLLARS

TRAUMA COLLARS LARGE



Large-sized grey Trauma collars.

Used to help prevent movement of the cervical spine and often are combined with lateral head blocks and straps.

[Read More](#)

SKU: CERV005

Price:

STIFF NECK COLLARS LARGE NO 4 NAVY BLUE



Stiff Neck Collars No 4 Navy Blue colour. Large size.


A neck collar helps support your head until the muscles and bones in your neck heal. Chronic neck pain has many causes — from things as minor as a backpack too full of textbooks to those as major as bone spurs or arthritis. In these cases, the collar helps relieve pain while treatment is ongoing.

[Read More](#)

SKU: CERV015

Price:

SOFT CERVICAL CREAM COLLARS LARGE



Large-sized cream soft cervical collars.

Used to help prevent movement of the cervical spine and often are combined with lateral head blocks and straps.

The soft Collar is made from comfortable foam to provide basic neck support. Lightweight, breathable and functional.

[Read More](#)

SKU: CERV001

Price:

INDEX

B

Bag only - Blue Technician Backpack 45x36x15cm (Eskom) 3

G

Gloves Latex High Risk Safety PF Large (50 pairs) 8

Gloves Latex PF Large (Box 100) 7

Gloves Latex Powdered Large (Box 100) 7

Gloves Latex Powdered Large (Pairs) 7

S

Safety pins Large Silver 45mm - 12s in bunch 5

Soft Cervical Cream Collars LARGE 12

Stiff Neck Collars LARGE No 4 Navy Blue 12

T

Tourniquet - Elastic & Buckle LARGE Black 10

Trauma Collars LARGE 12