

BJ TRADE+LINK



PRODUCT CATALOGUE 2023

TABLE OF CONTENTS

First Aid Kit Items	2
Spine Boards & Stretchers	4



FIRST AID KIT ITEMS



BLANKET FIRE/B IN POUCH 1 MTR X 1.2 MTR

Fire Blanket in Pouch 1 mtr x 1.2 mtr

[Read More](#)

SKU: BURN007

Price:



SPINE BOARDS & STRETCHERS

SPINE BOARD HEAVY-DUTY PLASTIC



Emergency Rescue Standard / Heavy-Duty Plastic (LDP) Spine board manufactured by BJ Tradelink.

Description: Lightweight Emergency Rescue Spine board with a unique tapered design to follow body contours ensuring easy manoeuvrability in confined spaces. Non-slip mild grain surface. Capable of floating in the water. Very strong and durable linear low-density UV Polyethylene plastic spine...

Emergency Rescue Standard / Heavy-Duty Plastic (LDP) Spine board manufactured by BJ Tradelink. Description: Lightweight Emergency Rescue Spine board with a unique tapered design to follow body contours ensuring easy maneuverability in confined spaces. Non-slip mild grain surface. Capable of floating in the water. Very strong and durable linear low-density UV Polyethylene plastic spine board. Impervious to all body fluids and easy to clean. Can be used in x-ray machines (100% X-ray translucent).

[Read More](#)

SKU: RESCU032

Price:



INDEX

B

Blanket Fire/B in Pouch 1 mtr x 1.2 mtr 3

S

Spine board Heavy-duty Plastic 5