

# BJ TRADE+LINK



## PRODUCT CATALOGUE 2023

# TABLE OF CONTENTS

**Consumable - Scissors/ Tweezers/ Forceps ..... 2**



## **CONSUMABLE - SCISSORS/ TWEEZERS/ FORCEPS**



## TWEEZER FORCEPS MOSQUITO STRAIGHT 12.5CM

TWEEZER Forceps Mosquito Straight 12.5cm.

A small, straight or curved forceps used in general surgery that has a locking grip with 3–5 teeth to allow ratchet clamping at various pressures. Mosquitos are used to retract in small fields, hold delicate tissue, and compress bleeding vessels, among other things

[Read More](#)

**SKU:** FCPS001

**Price:**



# INDEX

## T

TWEEZER Forceps Mosquito Straight 12.5cm 3