

# BJ TRADE+LINK



## PRODUCT CATALOGUE 2023

# TABLE OF CONTENTS

**Spine Boards & Stretchers** ..... 2



## **SPINE BOARDS & STRETCHERS**

## SPINE BOARD HEAVY-DUTY PLASTIC



Emergency Rescue Standard / Heavy-Duty Plastic (LDP) Spine board manufactured by BJ Tradelink.

Description: Lightweight Emergency Rescue Spine board with a unique tapered design to follow body contours ensuring easy manoeuvrability in confined spaces. Non-slip mild grain surface. Capable of floating in the water. Very strong and durable linear low-density UV Polyethylene plastic spine...

Emergency Rescue Standard / Heavy-Duty Plastic (LDP) Spine board manufactured by BJ Tradelink. Description: Lightweight Emergency Rescue Spine board with a unique tapered design to follow body contours ensuring easy maneuverability in confined spaces. Non-slip mild grain surface. Capable of floating in the water. Very strong and durable linear low-density UV Polyethylene plastic spine board. Impervious to all body fluids and easy to clean. Can be used in x-ray machines (100% X-ray translucent).

[Read More](#)

**SKU:** RESCU032

**Price:**





# INDEX

## S

Spine board Heavy-duty Plastic 3