

BJ TRADE+LINK



PRODUCT CATALOGUE 2023

TABLE OF CONTENTS

First Aid Bags	2
First Aid Kit Items	4
First Aid Kits	7
Spine boards / Stretcher Accessories	9
Spine Boards & Stretchers	12
Splints	15
Trauma Collars	17



FIRST AID BAGS



BAG ONLY - MOTORIST VEHICLE A4 BLACK

Bag only - A4 bag. Mostly used as Motor vehicle bags.

Colour: Black

Size: A4

[Read More](#)

SKU: BAGFA012

Price:



FIRST AID KIT ITEMS

RESUSAID CPR RED WITH CLEAR APRON



Resusaid Red CPR 165x165mm with clear apron.

A Resusaid or CPR mouthpiece, is a small device used to safely deliver rescue breaths during a cardiac arrest or respiratory arrest. Air is administered to the patient when the emergency responder exhales through the one-way filter valve. The CPR Mouthpiece has a built in one-way valve with a filter to protect the emergency responder from the...

[Read More](#)

SKU: RESC034

Price:



RESUSAID CPR 165X165MM WITH MILKY APRON



Resusaid CPR 165x165mm with milky apron.

A Resusaid or CPR mouthpiece, is a small device used to safely deliver rescue breaths during a cardiac arrest or respiratory arrest. Air is administered to the patient when the emergency responder exhales through the one-way filter valve. The CPR Mouthpiece has a built in one-way valve with a filter to protect the emergency responder from the...

[Read More](#)

SKU: RESC006

Price:

RESUSAID CPR 165X165MM WITH CLEAR APRON



Resusaid CPR 165x165mm with clear apron.

A Resusaid or CPR mouthpiece, is a small device used to safely deliver rescue breaths during a cardiac arrest or respiratory arrest. Air is administered to the patient when the emergency responder exhales through the one-way filter valve. The CPR Mouthpiece has a built in one-way valve with a filter to protect the emergency responder from the...

[Read More](#)

SKU: RESC023

Price:



FIRST AID KITS

Qty Description 1 x Resusaid 165x165mm Milky with instructions 1 x Gloves Nitrile Powder-Free (Large) pair 1 x Antiseptic Wipes (200mm x 190mm) [Read More](#)

SKU: BAGFA006

Price:

Qty Description 1 x Cetrimide wound cleaner 100ml 1 x Swabs Gauze 75x75x100 8ply n/s 2 x Swabs Gauze 75x75mm 5's Sterile 8ply 1 x Cotton Wool 100g N/Interleaved 1 x Tweezer S/Steel 85mm 1 x Scissors 16cm Bandage Black Handle+Lip 2 x Safety Pins ass 8's Bunch 2 x Pack of 8 safety pins (assorted sizes)(10cm) 2 x Gloves Latex PP Med n/s pairs 2 x Gloves Latex PP Large n/s pairs 2 x Splint Interlocking twisting green 4 x Conform Bandage i/w 75mmx 4.5m 12's 4 x Conform Bandage i/w 100mmx 4.5m 12's 1 x FAD No 3 Pad 75mm x 100mm x 2.25mt 1 x FAD No 5 Pad 200mmx150mm B100mmx2.5mtr 1 x Plaster Clear WP 100s Sterile 19x72mm 1 x Plaster Roll Elastic Beige 25mm x 3m 1 x Microporous Paper Tape 2.5cm x 5m 4 x Triangular Bandage N/W 40g 96x96x136 1 x Burn Gel Dressing 100x100 Sterile 1 x Eye Bath ribbed 'Plastic' [Read More](#)

SKU: FAKIT009

Price:

Qty Description 1 x Cetrimide wound cleaner 100ml 1 x Swabs Gauze 75x75x100 8ply n/s 2 x Swabs Gauze 75x75mm 5's Sterile 8ply 1 x Cotton Wool 100g N/Interleaved 1 x Tweezer S/Steel 85mm 1 x Scissors 16cm Bandage Black Handle+Lip 2 x Safety Pins ass 8's Bunch 2 x Pack of 8 safety pins (assorted sizes)(10cm) 2 x Gloves Latex PP Med n/s pairs 2 x Gloves Latex PP Large n/s pairs 2 x Splint Interlocking twisting green 4 x Conform Bandage i/w 75mmx 4.5m 12's 4 x Conform Bandage i/w 100mmx 4.5m 12's 1 x FAD No 3 Pad 75mm x 100mm x 2.25mt 1 x FAD No 5 Pad 200mmx150mm B100mmx2.5mtr 1 x Plaster Clear WP 100s Sterile 19x72mm 1 x Plaster Roll Elastic Beige 25mm x 3m 1 x Microporous Paper Tape 2.5cm x 5m 4 x Triangular Bandage N/W 40g 96x96x136 1 x Burn Gel Dressing 100x100 Sterile 1 x Eye Bath ribbed 'Plastic' [Read More](#)

SKU: BAGFA034

Price:

Qty Description 1 x Eyepad with single bandage n/s chett 1 x Cetrimide wound cleaner 50ml 1 x Swabs Gauze 75x75mm 5's Sterile 8ply 1 x Swabs gauze 75x75x8ply N/S pkt 10's 1 x Tweezer S/Steel 85mm 1 x Scissors 10cm pocket plastic handle red 2 x Safety Pins ass 12's Bunch 2 x Pack of 8 safety pins (assorted sizes)(10cm) 2 x Triangular Bandage N/W 40g 96x96x136 1 x Conform MT Lite 50mm x 4.5m Each 1 x Conform MT Lite 75mm x 4.5m Each 1 x Microporous Paper Tape 2.5cm x 5m 1 x Plaster Roll Elastic Beige 25mm x 3m 1 x Plaster Assorted 25's w/p/s beige 1 x FAD No 3 Pad 75mm x 100mm 1 x FAD No 5 Pad 200mmx150mm B100mmx2.5mtr 1 x Splint Combo disposable arm/leg 90deg 1 x Gloves Latex PP Medium n/s pairs 1 x Gloves Latex PP Large n/s pairs 1 x Resusaid 165x165mm f/m Quicksave yellow 1 x Blanket Emergency Rescue Adult 1 x Eye Bath ribbed 'Plastic' 1 x Cotton Wool 100gram Non-interleaved [Read More](#)

SKU: BAGFA001

Price:



SPINE BOARDS / STRETCHER ACCESSORIES

SPINE BOARD STRAPS WITH BUCKLES



Spine board or Stretcher Colour-coded adjustable disposable straps with quick-release buckles.

Used to secure patient on spine board or stretcher. Adjustable for different-sized persons. Quick release buckle for ease of securing the patient to the spine board. Colour-coded for ease of identifying the placement of straps around the spine board.

Colour-coded: Yellow, Orange, Black

3...

Spine board or Stretcher Colour-coded adjustable disposable straps with quick-release buckles. Used to secure patient on spine board or stretcher. Adjustable for different-sized persons. Quick release buckle for ease of securing the patient to the spine board. Colour-coded for ease of identifying the placement of straps around the spine board.

[Read More](#)

SKU: STRETCHER010

Price:



HEAD IMMOBILISER BOYD



Head Immobiliser Boyd

Easy to use and compatible with most spine boards on the market. The side blocks with wide Velcro closures, allow for a comfortable fit. It also comes with a head/chin straps and a base plate that fastens easily to the board. The Head Immobilizer has no metal components and can be worn during x-rays, MRI and CT scan procedures. Easy to clean and disinfect for...

Head Immobiliser Boyd. Easy to use and compatible with most spine boards on the market. The side blocks with wide Velcro closures, allow for a comfortable fit. It also comes with a head/chin straps and a base plate that fastens easily to the board. The Head Immobilizer has no metal components and can be worn during x-rays, MRI and CT scan procedures. Easy to clean and disinfect for reuse. The Head Immobilizer is made with the highest grade of foam.

[Read More](#)

SKU: RESC010

Price:



HEAD IMMOBILISER BLOCK ECONO

Trauma Head Immobiliser Block for adults "Econo"

Vinyl coated

Colour: Red

[Read More](#)

SKU: RESC026

Price:



HEAD IMMOBILISER BLOCK

Trauma Head Immobiliser Block for Adults.

Used on casualties with suspected head, neck or spinal injuries. Helps to prevent further injuries during movement and transportation. Easy to use, lightweight, the universal head immobiliser will not absorb bodily fluids, making it easy to clean and decontaminate between uses. Head supports are made of two plastic-coated, closed-cell foam with...

Trauma Head Immobiliser Block for Adults.

Used on casualties with suspected head, neck or spinal injuries. Helps to prevent further injuries during movement and transportation. Easy to use, lightweight, the universal head immobiliser will not absorb bodily fluids, making it easy to clean and decontaminate between uses. Head supports are made of two plastic-coated, closed-cell foam with large ear holes for monitoring the patient's ear canals. A universal attachment base and two durable, comfortable coated foam straps which fit snugly to the head and chin area. X-Ray/CT/MRI compatible and suitable for most extrication collars, cervical collars, spine boards and stretchers.

[Read More](#)

SKU: RESC003

Price:



SPINE BOARDS & STRETCHERS

STRETCHER PVC DOUBLE FOLD INCL. CARRY BAG



PVC Double Fold stretcher with aluminium legs including Carry Bag.

Aluminium alloy folding stretcher used for moving patients who require medical care. Easy storage. Lightweight and strong aluminium alloy frame. Integrated patient restraints to secure patients to the device while transporting.

Product dimensions (L×W×H): 220cm x 55cm x 15cm

Folded dimension: 94cm x 20cm x 9cm

Packing...

[Read More](#)

SKU: STRETCHER004

Price:



SPINE BOARD HEAVY-DUTY PLASTIC

Emergency Rescue Standard / Heavy-Duty Plastic (LDP) Spine board manufactured by BJ Tradelink.

Description: Lightweight Emergency Rescue Spine board with a unique tapered design to follow body contours ensuring easy manoeuvrability in confined spaces. Non-slip mild grain surface. Capable of floating in the water. Very strong and durable linear low-density UV Polyethylene plastic spine...

Emergency Rescue Standard / Heavy-Duty Plastic (LDP) Spine board manufactured by BJ Tradelink. Description: Lightweight Emergency Rescue Spine board with a unique tapered design to follow body contours ensuring easy maneuverability in confined spaces. Non-slip mild grain surface. Capable of floating in the water. Very strong and durable linear low-density UV Polyethylene plastic spine board. Impervious to all body fluids and easy to clean. Can be used in x-ray machines (100% X-ray translucent).



[Read More](#)

SKU: RESCU032

Price:



SPLINTS

SPLINT KRAMER (0.5M)

Kramer Wire Splint



The Kramer Wire Splint is designed to immobilise the injured body part and protects it from further injury. It aids in the reduction of pain, discomfort and swelling by limiting the movement of the injured limb. The Kramer Wire Splint is an extremely pliable multi-purpose splint and is built from a thin core of aluminium wire mesh, sandwiched between two layers of closed...

[Read More](#)

SKU: SP005

Price:



SPLINT COMBO ARM/LEG DISPOSABLE



Combo arm/leg 90-degree disposable splint.

An easily disposable and lightweight splint made of recyclable plastic. The combo splint can clip into a 180° and 90° configuration. Cost-effective solution. Can be used in first aid kits.

[Read More](#)

SKU: SPL002

Price:





TRAUMA COLLARS

TRAUMA SOFT AMBULANCE CERVICAL COLLARS WITH EARS (SMALL)



Trauma Soft Ambulance cervical collars with ears.

Size: Small.

Colour: Grey.

Used to help prevent movement of the cervical spine and often are combined with lateral head blocks and straps.

The soft Collar is made from comfortable foam to provide basic neck support. Lightweight, breathable and functional.

[Read More](#)

SKU: CERV018

Price:

TRAUMA COLLARS SMALL

Small grey Trauma collars.

Used to help prevent movement of the cervical spine and often are combined with lateral head blocks and straps.

[Read More](#)

SKU: CERV007

Price:

TRAUMA COLLARS MEDIUM

Medium-sized grey Trauma collars.

Used to help prevent movement of the cervical spine and often are combined with lateral head blocks and straps.

[Read More](#)

SKU: CERV006

Price:



TRAUMA COLLARS LARGE


Large-sized grey Trauma collars.

Used to help prevent movement of the cervical spine and often are combined with lateral head blocks and straps.

[Read More](#)

SKU: CERV005

Price:



SOFT CERVICAL CREAM COLLARS X-LARGE

Extra-Large size cream soft cervical collars.


Used to help prevent movement of the cervical spine and often are combined with lateral head blocks and straps.

The soft Collar is made from comfortable foam to provide basic neck support. Lightweight, breathable and functional.

[Read More](#)

SKU: CERV004

Price:



SOFT CERVICAL CREAM COLLARS SMALL

Small-sized cream soft cervical collars.


Used to help prevent movement of the cervical spine and often are combined with lateral head blocks and straps.

The soft Collar is made from comfortable foam to provide basic neck support. Lightweight, breathable and functional.

[Read More](#)

SKU: CERV003

Price:



SOFT CERVICAL CREAM COLLARS MEDIUM

Medium-sized cream soft cervical collars.


Used to help prevent movement of the cervical spine and often are combined with lateral head blocks and straps.

The soft Collar is made from comfortable foam to provide basic neck support. Lightweight, breathable and functional.

[Read More](#)

SKU: CERV002

Price:



SOFT CERVICAL CREAM COLLARS LARGE

Large-sized cream soft cervical collars.

Used to help prevent movement of the cervical spine and often are combined with lateral head blocks and straps.

The soft Collar is made from comfortable foam to provide basic neck support. Lightweight, breathable and functional.

[Read More](#)

SKU: CERV001

Price:

INDEX

B

Bag only - Motorist Vehicle A4 Black 3

H

Head Immobiliser Block 11

Head Immobiliser Block Econo 11

Head Immobiliser Boyd 10

R

Resusaid CPR 165x165mm with clear apron 6

Resusaid CPR 165x165mm with Milky apron 5

Resusaid CPR Red with Clear apron 5

S

Soft Cervical Cream Collars LARGE 20

Soft Cervical Cream Collars MEDIUM 20

Soft Cervical Cream Collars SMALL 20

Soft Cervical Cream Collars X-LARGE 19

Spine board Heavy-duty Plastic 13

Spine board Straps with buckles 10

Splint Combo arm/leg disposable 16

Splint Kramer (0.5m) 16

Stretcher PVC Double Fold incl. carry bag 13

T

Trauma Collars LARGE 19

Trauma Collars MEDIUM 18

Trauma Collars SMALL 18

Trauma Soft Ambulance Cervical Collars with ears (SMALL) 18