

# BJ TRADE+LINK



## PRODUCT CATALOGUE 2023

# TABLE OF CONTENTS

**Spine boards / Stretcher Accessories ..... 2**



## **SPINE BOARDS / STRETCHER ACCESSORIES**

## SPINE BOARD STRAPS WITH BUCKLES



Spine board or Stretcher Colour-coded adjustable disposable straps with quick-release buckles.

Used to secure patient on spine board or stretcher. Adjustable for different-sized persons. Quick release buckle for ease of securing the patient to the spine board. Colour-coded for ease of identifying the placement of straps around the spine board.

Colour-coded: Yellow, Orange, Black

3...

Spine board or Stretcher Colour-coded adjustable disposable straps with quick-release buckles. Used to secure patient on spine board or stretcher. Adjustable for different-sized persons. Quick release buckle for ease of securing the patient to the spine board. Colour-coded for ease of identifying the placement of straps around the spine board.

[Read More](#)

**SKU:** STRETCHER010

**Price:**



## SPIDER HARNESS TOP & TOE WITH VELCRO



Spider Harness Top & Toe with velcro fastening.

The Spider harness is an innovative strapping system for securing patients to the full spine board or scoop stretcher. The Spine board 12 Point Spider harness has 10 black straps for the torso and 2 red straps to easily identify to place at the top of the spine board for the upper body. These spider harness straps are made of heavy nylon...

Spider Harness Top & Toe with velcro fastening. The Spider harness is an innovative strapping system for securing patients to the full spine board or scoop stretcher. The Spine board 12 Point Spider harness has 10 black straps for the torso and 2 red straps to easily identify to place at the top of the spine board for the upper body. These spider harness straps are made of heavy nylon and are superior to other straps on the market. Black straps with red straps identifying the top of the spider harness.

[Read More](#)**SKU:** RESC001**Price:**

## SPIDER HARNESS TOP & TOE WITH BUCKLES

Spider Harness Top & Toe with buckle fastening.

The Spider harness is an innovative strapping system for securing patients to the full spine board or scoop stretcher. The Spine board 12 Point Spider harness has 10 black straps for the torso and 2 red straps to easily identify to place at the top of the spine board for the upper body. These spider harness straps are made of heavy nylon...

Spider Harness Top & Toe with buckles. Features buckle fastening. The Spider harness is an innovative strapping system for securing patients to the full spine board or scoop stretcher. The Spine board 12 Point Spider harness has 10 black straps for the torso and 2 red straps. This is to easily identify to place at the top of the spine board for the upper body. These spider harness straps are made of heavy nylon and are superior to other straps on the market. Black straps with red straps identifying the top of the spider harness.

[Read More](#)**SKU:** RESC016**Price:**

## SPIDER HARNESS COLOUR-CODED



Spider Harness Colour-coded Top and Toe strap.

A Spider Harness is used to immobilise a casualty with severe injuries while being transported on a spine board or stretcher. The harness one-piece is easily threaded through the spine board and secured around the patient's torso and limbs using high-strength Velcro.

This is the most reasonably priced colour-coded top-and-toe spider harness...

Spider Harness is Colour-coded Top and Toe strap. A Spider Harness is used to immobilise a casualty with severe injuries while being transported on a spine board or stretcher. The harness one-piece. It is easily threaded through the spine board and secured around the patient's torso and limbs using high-strength Velcro. This is the most reasonably priced colour-coded top-and-toe spider harness but does not compromise on quality.

[Read More](#)

**SKU:** RESC024

**Price:**



## SPIDER HARNESS BLACK TOP-ONLY



Spider Harness Black Top-only.

A Spider Harness is used to immobilise a casualty with severe injuries while being transported on a spine board or stretcher. The harness one-piece is easily threaded through the spine board and secured around the patient's torso and limbs using high-strength Velcro.

This is the most reasonably priced spider harness but does not compromise on...


Spider Harness Black Top-only. A Spider Harness is used to immobilize a casualty with severe injuries while being transported on a spine board or stretcher. The harness one-piece is easily threaded through the spine board and secured around the patient's torso and limbs using high-strength Velcro.

[Read More](#)

**SKU:** RESC025



**Price:**



## EXTRA-WIDTH SPIDER HARNESS TOP TO TOE WITH VELCRO

Spider Harness Top & Toe with Velcro fastening for extra-width patients.

The Spider harness is an innovative strapping system for securing patients to the full spine board or scoop stretcher. The Spine board 12 Point Spider harness is colour coded. These spider harness straps are made of heavy nylon and are superior to other straps on the market.

Extra-width for large to...

Spider Harness Top to Toe with Velcro fastening for extra-width patients. The Spider harness is an innovative strapping system for securing patients to the full spine board or scoop stretcher. The Spine board 12 Point Spider harness is colour coded. These spider harness straps are made of heavy nylon and are superior to other straps on the market. Extra-width for large to extra-large-sized patients.

[Read More](#)

**SKU:** RESCO041

**Price:**





# INDEX

## **E**

Extra-width Spider Harness Top to Toe with Velcro 6

## **S**

Spider Harness Black Top-only 5

Spider Harness Colour-Coded 5

Spider Harness Top & Toe with buckles 4

Spider Harness Top & Toe with Velcro 3

Spine board Straps with buckles 3