

BJ TRADE+LINK



PRODUCT CATALOGUE 2023

TABLE OF CONTENTS

Spine boards / Stretcher Accessories 2
Spine Boards & Stretchers 4



SPINE BOARDS / STRETCHER ACCESSORIES

SPINE BOARD STRAPS WITH BUCKLES



Spine board or Stretcher Colour-coded adjustable disposable straps with quick-release buckles.

Used to secure patient on spine board or stretcher. Adjustable for different-sized persons. Quick release buckle for ease of securing the patient to the spine board. Colour-coded for ease of identifying the placement of straps around the spine board.

Colour-coded: Yellow, Orange, Black

3...

Spine board or Stretcher Colour-coded adjustable disposable straps with quick-release buckles. Used to secure patient on spine board or stretcher. Adjustable for different-sized persons. Quick release buckle for ease of securing the patient to the spine board. Colour-coded for ease of identifying the placement of straps around the spine board.

[Read More](#)

SKU: STRETCHER010

Price:





SPINE BOARDS & STRETCHERS

STRETCHER PVC DOUBLE FOLD INCL. CARRY BAG



PVC Double Fold stretcher with aluminium legs including Carry Bag.

Aluminium alloy folding stretcher used for moving patients who require medical care. Easy storage. Lightweight and strong aluminium alloy frame. Integrated patient restraints to secure patients to the device while transporting.

Product dimensions (L×W×H): 220cm x 55cm x 15cm

Folded dimension: 94cm x 20cm x 9cm

Packing...

[Read More](#)

SKU: STRETCHER004

Price:



SPINE BOARD HEAVY-DUTY PLASTIC

Emergency Rescue Standard / Heavy-Duty Plastic (LDP) Spine board manufactured by BJ Tradelink.

Description: Lightweight Emergency Rescue Spine board with a unique tapered design to follow body contours ensuring easy manoeuvrability in confined spaces. Non-slip mild grain surface. Capable of floating in the water. Very strong and durable linear low-density UV Polyethylene plastic spine...

Emergency Rescue Standard / Heavy-Duty Plastic (LDP) Spine board manufactured by BJ Tradelink. Description: Lightweight Emergency Rescue Spine board with a unique tapered design to follow body contours ensuring easy maneuverability in confined spaces. Non-slip mild grain surface. Capable of floating in the water. Very strong and durable linear low-density UV Polyethylene plastic spine board. Impervious to all body fluids and easy to clean. Can be used in x-ray machines (100% X-ray translucent).



[Read More](#)

SKU: RESCU032

Price:

INDEX

S

Spine board Heavy-duty Plastic 5
Spine board Straps with buckles 3
Stretcher PVC Double Fold incl. carry bag 5