

BJ TRADE+LINK



PRODUCT CATALOGUE 2023

TABLE OF CONTENTS

First Aid Bags	2
First Aid Boxes	4
Spine boards / Stretcher Accessories	6
Spine Boards & Stretchers	11
Splints	13
Trauma Collars	15



FIRST AID BAGS



BAG ONLY - MOTORIST VEHICLE A4 BLACK

Bag only - A4 bag. Mostly used as Motor vehicle bags.

Colour: Black

Size: A4

[Read More](#)

SKU: BAGFA012

Price:



FIRST AID BOXES

FIRST AID BOX LDP GREEN LARGE



Large Plastic Green First Aid Box made with UV Polyethylene material manufactured by BJ Tradelink, the original supplier.

Very durable bottle green plastic large box used mainly for first aid kits in the mines and construction sites with 6 compartments secured by metal dividers. Straps on the outside for fleece blankets.

Material: UV Polyethylene

Wall Thickness: ...

[Read More](#)

SKU: LDP001

Price:





SPINE BOARDS / STRETCHER ACCESSORIES

SPINE BOARD STRAPS WITH BUCKLES



Spine board or Stretcher Colour-coded adjustable disposable straps with quick-release buckles.

Used to secure patient on spine board or stretcher. Adjustable for different-sized persons. Quick release buckle for ease of securing the patient to the spine board. Colour-coded for ease of identifying the placement of straps around the spine board.

Colour-coded: Yellow, Orange, Black

3...

Spine board or Stretcher Colour-coded adjustable disposable straps with quick-release buckles. Used to secure patient on spine board or stretcher. Adjustable for different-sized persons. Quick release buckle for ease of securing the patient to the spine board. Colour-coded for ease of identifying the placement of straps around the spine board.

[Read More](#)

SKU: STRETCHER010

Price:



SPIDER HARNESS TOP & TOE WITH VELCRO



Spider Harness Top & Toe with velcro fastening.

The Spider harness is an innovative strapping system for securing patients to the full spine board or scoop stretcher. The Spine board 12 Point Spider harness has 10 black straps for the torso and 2 red straps to easily identify to place at the top of the spine board for the upper body. These spider harness straps are made of heavy nylon...

Spider Harness Top & Toe with velcro fastening. The Spider harness is an innovative strapping system for securing patients to the full spine board or scoop stretcher. The Spine board 12 Point Spider harness has 10 black straps for the torso and 2 red straps to easily identify to place at the top of the spine board for the upper body. These spider harness straps are made of heavy nylon and are superior to other straps on the market. Black straps with red straps identifying the top of the spider harness.


[Read More](#)**SKU:** RESC001**Price:**

SPIDER HARNESS TOP & TOE WITH BUCKLES

Spider Harness Top & Toe with buckle fastening.

The Spider harness is an innovative strapping system for securing patients to the full spine board or scoop stretcher. The Spine board 12 Point Spider harness has 10 black straps for the torso and 2 red straps to easily identify to place at the top of the spine board for the upper body. These spider harness straps are made of heavy nylon...

Spider Harness Top & Toe with buckles. Features buckle fastening. The Spider harness is an innovative strapping system for securing patients to the full spine board or scoop stretcher. The Spine board 12 Point Spider harness has 10 black straps for the torso and 2 red straps. This is to easily identify to place at the top of the spine board for the upper body. These spider harness straps are made of heavy nylon and are superior to other straps on the market. Black straps with red straps identifying the top of the spider harness.

[Read More](#)**SKU:** RESC016**Price:**

HEAD IMMOBILISER BOYD



Head Immobiliser Boyd

Easy to use and compatible with most spine boards on the market. The side blocks with wide Velcro closures, allow for a comfortable fit. It also comes with a head/chin straps and a base plate that fastens easily to the board. The Head Immobilizer has no metal components and can be worn during x-rays, MRI and CT scan procedures. Easy to clean and disinfect for...

Head Immobiliser Boyd. Easy to use and compatible with most spine boards on the market. The side blocks with wide Velcro closures, allow for a comfortable fit. It also comes with a head/chin straps and a base plate that fastens easily to the board. The Head Immobilizer has no metal components and can be worn during x-rays, MRI and CT scan procedures. Easy to clean and disinfect for reuse. The Head Immobilizer is made with the highest grade of foam.

[Read More](#)

SKU: RESC010

Price:

HEAD IMMOBILISER BLOCK ECONO



Trauma Head Immobiliser Block for adults "Econo"

Vinyl coated

Colour: Red

[Read More](#)

SKU: RESC026

Price:

HEAD IMMOBILISER BLOCK

Trauma Head Immobiliser Block for Adults.

Used on casualties with suspected head, neck or spinal injuries. Helps to prevent further injuries during movement and transportation. Easy to use, lightweight, the universal head immobiliser will not absorb bodily fluids,

making it easy to clean and decontaminate between uses. Head supports are made of two plastic-coated, closed-cell foam with...



Trauma Head Immobiliser Block for Adults.

Used on casualties with suspected head, neck or spinal injuries. Helps to prevent further injuries during movement and transportation. Easy to use, lightweight, the universal head immobiliser will not absorb bodily fluids, making it easy to clean and decontaminate between uses. Head supports are made of two plastic-coated, closed-cell foam with large ear holes for monitoring the patient's ear canals. A universal attachment base and two durable, comfortable coated foam straps which fit snugly to the head and chin area. X-Ray/CT/MRI compatible and suitable for most extrication collars, cervical collars, spine boards and stretchers.

[Read More](#)

SKU: RESC003

Price:

EXTRA-WIDTH SPIDER HARNESS TOP TO TOE WITH VELCRO



Spider Harness Top & Toe with Velcro fastening for extra-width patients.

The Spider harness is an innovative strapping system for securing patients to the full spine board or scoop stretcher. The Spine board 12 Point Spider harness is colour coded. These spider harness straps are made of heavy nylon and are superior to other straps on the market.

Extra-width for large to...

Spider Harness Top to Toe with Velcro fastening for extra-width patients. The Spider harness is an innovative strapping system for securing patients to the full spine board or scoop stretcher. The Spine board 12 Point Spider harness is colour coded. These spider harness straps are made of heavy nylon and are superior to other straps on the market. Extra-width for large to extra-large-sized patients.

[Read More](#)

SKU: RESC0041

Price:





SPINE BOARDS & STRETCHERS

SPINE BOARD HEAVY-DUTY PLASTIC



Emergency Rescue Standard / Heavy-Duty Plastic (LDP) Spine board manufactured by BJ Tradelink.

Description: Lightweight Emergency Rescue Spine board with a unique tapered design to follow body contours ensuring easy manoeuvrability in confined spaces. Non-slip mild grain surface. Capable of floating in the water. Very strong and durable linear low-density UV Polyethylene plastic spine...

Emergency Rescue Standard / Heavy-Duty Plastic (LDP) Spine board manufactured by BJ Tradelink. Description: Lightweight Emergency Rescue Spine board with a unique tapered design to follow body contours ensuring easy maneuverability in confined spaces. Non-slip mild grain surface. Capable of floating in the water. Very strong and durable linear low-density UV Polyethylene plastic spine board. Impervious to all body fluids and easy to clean. Can be used in x-ray machines (100% X-ray translucent).

[Read More](#)

SKU: RESCU032

Price:





SPLINTS

SPLINT KRAMER (0.5M)

Kramer Wire Splint



The Kramer Wire Splint is designed to immobilise the injured body part and protects it from further injury. It aids in the reduction of pain, discomfort and swelling by limiting the movement of the injured limb. The Kramer Wire Splint is an extremely pliable multi-purpose splint and is built from a thin core of aluminium wire mesh, sandwiched between two layers of closed...

[Read More](#)

SKU: SP005

Price:



SPLINT COMBO ARM/LEG DISPOSABLE



Combo arm/leg 90-degree disposable splint.

An easily disposable and lightweight splint made of recyclable plastic. The combo splint can clip into a 180° and 90° configuration. Cost-effective solution. Can be used in first aid kits.

[Read More](#)

SKU: SPL002

Price:





TRAUMA COLLARS

TRAUMA SOFT AMBULANCE CERVICAL COLLARS WITH EARS (SMALL)



Trauma Soft Ambulance cervical collars with ears.

Size: Small.

Colour: Grey.

Used to help prevent movement of the cervical spine and often are combined with lateral head blocks and straps.

The soft Collar is made from comfortable foam to provide basic neck support. Lightweight, breathable and functional.

[Read More](#)

SKU: CERV018

Price:

TRAUMA COLLARS SMALL

Small grey Trauma collars.

Used to help prevent movement of the cervical spine and often are combined with lateral head blocks and straps.

[Read More](#)

SKU: CERV007

Price:

TRAUMA COLLARS MEDIUM

Medium-sized grey Trauma collars.

Used to help prevent movement of the cervical spine and often are combined with lateral head blocks and straps.

[Read More](#)

SKU: CERV006

Price:



TRAUMA COLLARS LARGE

Large-sized grey Trauma collars.

Used to help prevent movement of the cervical spine and often are combined with lateral head blocks and straps.

[Read More](#)

SKU: CERV005

Price:



SOFT CERVICAL CREAM COLLARS X-LARGE

Extra-Large size cream soft cervical collars.

Used to help prevent movement of the cervical spine and often are combined with lateral head blocks and straps.

The soft Collar is made from comfortable foam to provide basic neck support. Lightweight, breathable and functional.

[Read More](#)

SKU: CERV004

Price:



SOFT CERVICAL CREAM COLLARS SMALL

Small-sized cream soft cervical collars.

Used to help prevent movement of the cervical spine and often are combined with lateral head blocks and straps.

The soft Collar is made from comfortable foam to provide basic neck support. Lightweight, breathable and functional.

[Read More](#)

SKU: CERV003

Price:



SOFT CERVICAL CREAM COLLARS MEDIUM

Medium-sized cream soft cervical collars.

Used to help prevent movement of the cervical spine and often are combined with lateral head blocks and straps.

The soft Collar is made from comfortable foam to provide basic neck support. Lightweight, breathable and functional.

[Read More](#)

SKU: CERV002

Price:



SOFT CERVICAL CREAM COLLARS LARGE

Large-sized cream soft cervical collars.

Used to help prevent movement of the cervical spine and often are combined with lateral head blocks and straps.

The soft Collar is made from comfortable foam to provide basic neck support. Lightweight, breathable and functional.

[Read More](#)

SKU: CERV001

Price:

INDEX

B

Bag only - Motorist Vehicle A4 Black 3

E

Extra-width Spider Harness Top to Toe with Velcro 10

F

First Aid Box LDP Green Large 5

H

Head Immobiliser Block 9

Head Immobiliser Block Econo 9

Head Immobiliser Boyd 9

S

Soft Cervical Cream Collars LARGE 18

Soft Cervical Cream Collars MEDIUM 18

Soft Cervical Cream Collars SMALL 18

Soft Cervical Cream Collars X-LARGE 17

Spider Harness Top & Toe with buckles 8

Spider Harness Top & Toe with Velcro 7

Spine board Heavy-duty Plastic 12

Spine board Straps with buckles 7

Splint Combo arm/leg disposable 14

Splint Kramer (0.5m) 14

T

Trauma Collars LARGE 17

Trauma Collars MEDIUM 16

Trauma Collars SMALL 16

Trauma Soft Ambulance Cervical Collars with ears (SMALL) 16