

# BJ TRADE+LINK



## PRODUCT CATALOGUE 2023

# TABLE OF CONTENTS

**Consumable - Scissors/ Tweezers/ Forceps ..... 2**



## **CONSUMABLE - SCISSORS/ TWEEZERS/ FORCEPS**

## TWEEZER PLASTIC



TWEEZER plastic is used in first aid kits or for various other uses.

[Read More](#)

**SKU:** TWEEZER002

**Price:**



## TWEEZER FORCEPS MOSQUITO STRAIGHT 12.5CM

TWEEZER Forceps Mosquito Straight 12.5cm.

A small, straight or curved forceps used in general surgery that has a locking grip with 3-5 teeth to allow ratchet clamping at various pressures. Mosquitos are used to retract in small fields, hold delicate tissue, and compress bleeding vessels, among other things

[Read More](#)

**SKU:** FCPS001

**Price:**

## TWEEZER FLAT STAINLESS STEEL METAL 88MM



TWEEZER Flat Stainless Steel Metal 88mm

[Read More](#)

**SKU:** TWEEZER001

**Price:**





# INDEX

## T

TWEEZER Flat Stainless Steel Metal 88mm 3

TWEEZER Forceps Mosquito Straight 12.5cm 3

TWEEZER Plastic 3