

BJ TRADE+LINK



PRODUCT CATALOGUE 2023

TABLE OF CONTENTS

Tourniquets 2



TOURNIQUETS



TOURNIQUET - TWISTING STICK 200MM X 10MM

Tourniquet - Tourniquet - Twisting stick 200mm x 10mm

Used in mining first aid kits.

[Read More](#)

SKU: SPL013

Price:



TOURNIQUET - TWISTING STICK 150MM X 10MM

Tourniquet - Tourniquet - Twisting stick 150mm x 10mm

Used in mining first aid kits.

[Read More](#)

SKU: SPL005

Price:



TOURNIQUET - QUICK RELEASE CLIP ROYAL BLUE



Tourniquet - Quick Release clip Royal Blue

The tourniquet allows for pressure to be applied to the arm so that venous blood returning to the heart can be slowed down. As a result, the blood vessel walls become temporarily occluded and the veins distend due to the pooling of blood. This allows veins to become more visible and easier to palpate.

How to use:

1.Quick release...

[Read More](#)

SKU: TOURN004

Price:

INDEX

T

Tourniquet - Quick Release clip Royal Blue 4

Tourniquet - Twisting stick 150mm x 10mm 3

Tourniquet - Twisting stick 200mm x 10mm 3